



Proteogenix

Proteogenix SAS

www.proteogenix.science
19 rue de la Haye
67300 Schiltigheim, France
Tel +33 (0) 390 20 54 70
Fax +33 (0) 978 53 36 90
products@proteogenix.fr

Anti- α -Smooth Muscle Actin Mouse Monoclonal Antibody (6A12)

OVERVIEW

PRODUCT REFERENCE	PTX1038
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1

PRODUCT INFORMATION

EXPRESSION SYSTEM	Mouse
BUFFER	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide
STORAGE CONDITIONS	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40% glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.



WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.