



Anti-GAPDH Rabbit Polyclonal Antibody(Zebrafish Specific)

OVERVIEW

PRODUCT REFERENCE	PTX1054
CLONALITY	Polyclonal Antibody
ISOTYPE	IgG

PRODUCT INFORMATION

EXPRESSION SYSTEM	Rabbit
BUFFER	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide
STORAGE CONDITIONS	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40% glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.

