



Anti- β II tubulin Mouse Monoclonal Antibody (Rodent Preferred)

OVERVIEW

PRODUCT REFERENCE	PTX1100
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1

PRODUCT INFORMATION

EXPRESSION SYSTEM	Mouse
BUFFER	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide
STORAGE CONDITIONS	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40% glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.

