



FITC* Goat Anti-Rabbit IgG (H+L)

OVERVIEW

PRODUCT REFERENCE	PTX17819
CLONALITY	Polyclonal Antibody
TARGET	Anti-rabbit antibodies

PRODUCT INFORMATION

HOST	Goat
REACTIVITY	Rabbit
BUFFER & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

