



# Proteogenix

**Proteogenix SAS**

www.proteogenix.science  
19 rue de la Haye  
67300 Schiltigheim, France

Tel +33 (0) 390 20 54 70  
Fax +33 (0) 978 53 36 90  
products@proteogenix.fr

## BIOTIN\* Goat Anti-Rabbit IgG (H+L)

### OVERVIEW

---

<b>PRODUCT REFERENCE</b>	PTX17821
<b>CLONALITY</b>	Polyclonal Antibody
<b>TARGET</b>	Anti-rabbit antibodies

### PRODUCT INFORMATION

---

<b>HOST</b>	Goat
<b>REACTIVITY</b>	Rabbit
<b>BUFFER &amp; STORAGE</b>	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles



**WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.**