



## HRP-conjugated anti-M13 phage antibody

### OVERVIEW

---

<b>PRODUCT REFERENCE</b>	PTX17839
<b>CLONALITY</b>	Monoclonal Antibody
<b>ISOTYPE</b>	IgG1
<b>TARGET</b>	M13 Bacteriophage

### PRODUCT INFORMATION

---

<b>HOST</b>	Mouse
<b>REACTIVITY</b>	M13 Bacteriophage
<b>FORM &amp; BUFFER</b>	Supplied as solution form in PBS, pH7.4, containing 0.02% NaN <sub>3</sub> , 50% glycerol.
<b>STORAGE</b>	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.

