



Anti MBP tag mouse monoclonal antibody

OVERVIEW

PRODUCT REFERENCE	PTX17862
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1

PRODUCT INFORMATION

EXPRESSION SYSTEM	Mouse
BUFFER	PBS pH7.4, containing 0.02% NaN ₃ and 50% glycerol
STORAGE CONDITIONS	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40% glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.

