



Proteogenix

Proteogenix SAS

www.proteogenix.science
19 rue de la Haye
67300 Schiltigheim, France

Tel +33 (0) 390 20 54 70
Fax +33 (0) 978 53 36 90
products@proteogenix.fr

Anti HA tag mouse monoclonal antibody(HRP)

OVERVIEW

| | |
|--------------------------|---------------------|
| PRODUCT REFERENCE | PTX17870 |
| CLONALITY | Monoclonal Antibody |
| ISOTYPE | IgG1 |

PRODUCT INFORMATION

| | |
|---------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| EXPRESSION SYSTEM | Mouse |
| BUFFER | PBS pH7.4, containing 0.02% NaN ₃ and 50% glycerol |
| STORAGE CONDITIONS | Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40% glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles. |



WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.