



Anti Cytokeratin 19 Monoclonal Antibody

OVERVIEW

PRODUCT REFERENCE	PTX17886
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1
TARGET	Synthesized CK19.

PRODUCT INFORMATION

HOST	Mouse
REACTIVITY	human
FORM & BUFFER	PBS, pH7.4, containing 0.05% proclin300, 50% glycerol.
STORAGE	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.

