



Proteogenix

Proteogenix SAS

www.proteogenix.science
19 rue de la Haye
67300 Schiltigheim, France

Tel +33 (0) 390 20 54 70
Fax +33 (0) 978 53 36 90
products@proteogenix.fr

Human IgG4 (S228P) Isotype Control antibody (HyHEL-10)

OVERVIEW

PRODUCT REFERENCE	PTX17904
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG4 (S228P)

PRODUCT INFORMATION

EXPRESSION SYSTEM	Human
BUFFER	PBS pH7.5
STORAGE CONDITIONS	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40% glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.



WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.