



## CD274 / PD-L1 / B7-H1 Monoclonal Monoclonal Antibody 2C11

### OVERVIEW

---

<b>PRODUCT REFERENCE</b>	PTX18518
<b>CLONALITY</b>	Monoclonal Antibody
<b>ISOTYPE</b>	IgG1
<b>TARGET</b>	CD274

### PRODUCT INFORMATION

---

<b>HOST</b>	Mouse
<b>REACTIVITY</b>	Human
<b>FORM &amp; BUFFER</b>	Liquid in PBS buffer, pH 7.4
<b>STORAGE</b>	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.

