



CD274 / PD-L1 / B7-H1 Monoclonal Monoclonal Antibody 5H8

OVERVIEW

PRODUCT REFERENCE	PTX18519
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG2b
TARGET	CD274

PRODUCT INFORMATION

HOST	Mouse
REACTIVITY	Human
FORM & BUFFER	Liquid in PBS buffer, pH 7.4
STORAGE	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.

