

## CXCL9 / MIG Polyclonal Polyclonal Antibody

### OVERVIEW

---

PRODUCT REFERENCE	PTX18552
CLONALITY	Polyclonal Antibody
ISOTYPE	IgG
TARGET	CXCL9

### PRODUCT INFORMATION

---

HOST	Rabbit
REACTIVITY	Human
FORM & BUFFER	Liquid in PBS buffer, pH 7.4
STORAGE	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.

