



HCV NS1/gp68/gp70/Envelope glycoprotein E2 Monoclonal Antibody HC33.4

OVERVIEW

PRODUCT REFERENCE	PTX18581
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1-kappa
TARGET	HCV NS1

PRODUCT INFORMATION

HOST	Mouse
REACTIVITY	Hepatitis C virus genotype 1a (isolate H77) (HCV)
FORM & BUFFER	Liquid in PBS buffer, pH 7.4
STORAGE	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.

