



# Proteogenix

**Proteogenix SAS**

www.proteogenix.science  
19 rue de la Haye  
67300 Schiltigheim, France  
Tel +33 (0) 390 20 54 70  
Fax +33 (0) 978 53 36 90  
products@proteogenix.fr

## HCV NS1/gp68/gp70/Envelope glycoprotein E2 Monoclonal Antibody HC84-27

### OVERVIEW

---

<b>PRODUCT REFERENCE</b>	PTX18585
<b>CLONALITY</b>	Monoclonal Antibody
<b>ISOTYPE</b>	IgG1-kappa
<b>TARGET</b>	HCV NS1

### PRODUCT INFORMATION

---

<b>HOST</b>	Mouse
<b>REACTIVITY</b>	Hepatitis C virus genotype 1a (isolate H77) (HCV)
<b>FORM &amp; BUFFER</b>	Liquid in PBS buffer, pH 7.4
<b>STORAGE</b>	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.



**WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.**