



## LAMTOR3 Polyclonal Polyclonal Antibody

### OVERVIEW

---

|                          |                     |
|--------------------------|---------------------|
| <b>PRODUCT REFERENCE</b> | PTX18623            |
| <b>CLONALITY</b>         | Polyclonal Antibody |
| <b>ISOTYPE</b>           | IgG                 |
| <b>TARGET</b>            | LAMTOR3             |

### PRODUCT INFORMATION

---

|                          |   |
|--------------------------|---|
| <b>HOST</b>              | Rabbit  |
| <b>REACTIVITY</b>        | Human   |
| <b>FORM &amp; BUFFER</b> | Liquid in PBS buffer, pH 7.4  |
| <b>STORAGE</b>           | Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles. |

