



SERPINB1 Polyclonal Polyclonal Antibody

OVERVIEW

| | |
|--------------------------|---------------------|
| PRODUCT REFERENCE | PTX18671 |
| CLONALITY | Polyclonal Antibody |
| ISOTYPE | IgG |
| TARGET | SERPINB1 |

PRODUCT INFORMATION

| | |
|--------------------------|---|
| HOST | Rabbit |
| REACTIVITY | Human |
| FORM & BUFFER | Liquid in PBS buffer, pH 7.4 |
| STORAGE | Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles. |

