



Proteogenix

Proteogenix SAS

www.proteogenix.science
19 rue de la Haye
67300 Schiltigheim, France
Tel +33 (0) 390 20 54 70
Fax +33 (0) 978 53 36 90
products@proteogenix.fr

TNF α / TNF-alpha Monoclonal Antibody CDP571

OVERVIEW

PRODUCT REFERENCE	PTX18819
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG4-kappa
TARGET	TNF α

PRODUCT INFORMATION

HOST	Humanized
REACTIVITY	Human
FORM & BUFFER	Liquid in PBS buffer, pH 7.4
STORAGE	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.



WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.