



Proteogenix

Proteogenix SAS

www.proteogenix.science
19 rue de la Haye
67300 Schiltigheim, France
Tel +33 (0) 390 20 54 70
Fax +33 (0) 978 53 36 90
products@proteogenix.fr

Monkeypox virus/MPXV M1R SAA0283 Monoclonal Antibody SAA0283

OVERVIEW

PRODUCT REFERENCE	PTX18821
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG2a, kappa
TARGET	M1R

PRODUCT INFORMATION

HOST	Mouse
REACTIVITY	Monkeypox virus(MPXV)
FORM & BUFFER	Liquid in PBS buffer, pH 7.4
STORAGE	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.



WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.