



Monkeypox virus/MPXV A29L Monoclonal Antibody SAA0286

OVERVIEW

PRODUCT REFERENCE	PTX18824
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG2a, lambda
TARGET	A29L

PRODUCT INFORMATION

HOST	Mouse
REACTIVITY	Monkeypox virus(MPXV)
FORM & BUFFER	Liquid in PBS buffer, pH 7.4
STORAGE	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.

