



Anti-Vaccinia H3L/Ag35 Antibody (B9)

OVERVIEW

PRODUCT REFERENCE	PTX18826
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG2a, kappa
TARGET	H3L/Ag35

PRODUCT INFORMATION

HOST	Mouse
REACTIVITY	Vaccinia
FORM & BUFFER	Liquid in PBS buffer, pH 7.4
STORAGE	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.

