



## Anti-Vaccinia A27 Antibody (8E3#)

### OVERVIEW

---

<b>PRODUCT REFERENCE</b>	PTX18830
<b>CLONALITY</b>	Monoclonal Antibody
<b>ISOTYPE</b>	IgG2a, lambda
<b>TARGET</b>	A27

### PRODUCT INFORMATION

---

<b>HOST</b>	Mouse
<b>REACTIVITY</b>	Vaccinia
<b>FORM &amp; BUFFER</b>	Liquid in PBS buffer, pH 7.4
<b>STORAGE</b>	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.

