



scFv isotype control antibody (HyHEL-10)

OVERVIEW

PRODUCT REFERENCE	PTX18838
CLONALITY	Monoclonal Antibody
ISOTYPE	scFv

PRODUCT INFORMATION

EXPRESSION SYSTEM

BUFFER PBS pH7.4, containing 0.02% NaN₃ and 50% glycerol

STORAGE CONDITIONS Store at 4 °C for short term. Freeze and store at -80 °C for long term.
Addition of 20 to 40% glycerol may improve cryoprotection.
Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions.
Aliquot the product according to your needs to avoid freeze/thaw cycles.

