



Proteogenix

Proteogenix SAS

www.proteogenix.science
19 rue de la Haye
67300 Schiltigheim, France
Tel +33 (0) 390 20 54 70
Fax +33 (0) 978 53 36 90
products@proteogenix.fr

Anti Strep tag mouse monoclonal antibody (C23.21)

OVERVIEW

PRODUCT REFERENCE	PTX18839
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1, kappa

PRODUCT INFORMATION

EXPRESSION SYSTEM	Mouse
BUFFER	PBS pH7.4, containing 0.02% NaN ₃ and 50% glycerol.
STORAGE CONDITIONS	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40% glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.



WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.