



Proteogenix

Proteogenix SAS

www.proteogenix.science
19 rue de la Haye
67300 Schiltigheim, France

Tel +33 (0) 390 20 54 70
Fax +33 (0) 978 53 36 90
products@proteogenix.fr

Anti-Monkeypox virus/MPXV A44R Polyclonal Antibody

OVERVIEW

PRODUCT REFERENCE	PTX18847
CLONALITY	Polyclonal Antibody
ISOTYPE	IgG
TARGET	A44R

PRODUCT INFORMATION

HOST	Rabbit
REACTIVITY	Monkeypox virus(MPXV)
FORM & BUFFER	Liquid in PBS buffer, pH 7.4
STORAGE	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.



WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.