



Proteogenix

Proteogenix SAS

www.proteogenix.science
19 rue de la Haye
67300 Schiltigheim, France

Tel +33 (0) 390 20 54 70
Fax +33 (0) 978 53 36 90
products@proteogenix.fr

Anti-Monkeypox virus/MPXV C18L Polyclonal Antibody

OVERVIEW

| | |
|--------------------------|---------------------|
| PRODUCT REFERENCE | PTX18848 |
| CLONALITY | Polyclonal Antibody |
| ISOTYPE | IgG |
| TARGET | C18L |

PRODUCT INFORMATION

| | |
|--------------------------|---|
| HOST | Rabbit |
| REACTIVITY | Monkeypox virus(MPXV) |
| FORM & BUFFER | Liquid in PBS buffer, pH 7.4 |
| STORAGE | Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles. |



WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.