



Anti-Monkeypox virus/MPXV D14L Polyclonal Antibody

OVERVIEW

PRODUCT REFERENCE	PTX18851
CLONALITY	Polyclonal Antibody
ISOTYPE	IgG
TARGET	D14L

PRODUCT INFORMATION

HOST	Rabbit
REACTIVITY	Monkeypox virus(MPXV)
FORM & BUFFER	Liquid in PBS buffer, pH 7.4
STORAGE	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.

