



Proteogenix

Proteogenix SAS

www.proteogenix.science
19 rue de la Haye
67300 Schiltigheim, France
Tel +33 (0) 390 20 54 70
Fax +33 (0) 978 53 36 90
products@proteogenix.fr

Anti-Denosumab Polyclonal Antibody

OVERVIEW

| | |
|--------------------------|---------------------|
| PRODUCT REFERENCE | PTX18888 |
| CLONALITY | Polyclonal Antibody |
| ISOTYPE | IgG |
| TARGET | Denosumab |

PRODUCT INFORMATION

| | |
|--------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------|
| HOST | Rabbit |
| REACTIVITY | Drug |
| FORM & BUFFER | 0.01M PBS, pH 7.4, 50% glycerol, 0.09% sodium azide. |
| STORAGE | Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles. |



WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.