



Anti-Cetuximab Polyclonal Antibody

OVERVIEW

PRODUCT REFERENCE	PTX18889
CLONALITY	Polyclonal Antibody
ISOTYPE	IgG
TARGET	Cetuximab

PRODUCT INFORMATION

HOST	Rabbit
REACTIVITY	Drug
FORM & BUFFER	0.01M PBS, pH 7.4, 50% glycerol, 0.09% sodium azide.
STORAGE	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.

