



Anti-Human EGFR/ERBB1/HER1 Nanobody (SAA1038)

OVERVIEW

PRODUCT REFERENCE	PTX18906
CLONALITY	Monoclonal Antibody
ISOTYPE	VHH-8His-Cys-tag
TARGET	EGFR/ERBB1/HER1

PRODUCT INFORMATION

HOST	XtenCHO
REACTIVITY	Human
FORM & BUFFER	0.01M PBS, pH 7.4.
STORAGE	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.

