



Anti-Human C4A Nanobody (SAA0908)

OVERVIEW

PRODUCT REFERENCE	PTX18934
CLONALITY	Monoclonal Antibody
ISOTYPE	VHH-8His-Cys-tag
TARGET	C4A

PRODUCT INFORMATION

HOST	XtenCHO
REACTIVITY	Human
FORM & BUFFER	0.01M PBS, pH 7.4.
STORAGE	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.

