



Anti-Human IL23A/IL-23p19 Nanobody (SAA1160)

OVERVIEW

PRODUCT REFERENCE	PTX18986
CLONALITY	Monoclonal Antibody
ISOTYPE	VHH-8His-Cys-tag
TARGET	IL23A/IL-23p19

PRODUCT INFORMATION

HOST	XtenCHO
REACTIVITY	Human
FORM & BUFFER	0.01M PBS, pH 7.4.
STORAGE	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.

