



## Anti-SARS-CoV-2 S Protein Nanobody (SAA0886)

### OVERVIEW

---

|                          |                     |
|--------------------------|---------------------|
| <b>PRODUCT REFERENCE</b> | PTX18994            |
| <b>CLONALITY</b>         | Monoclonal Antibody |
| <b>ISOTYPE</b>           | VHH-8His-Cys-tag    |
| <b>TARGET</b>            | S Protein           |

### PRODUCT INFORMATION

---

|                          |   |
|--------------------------|---|
| <b>HOST</b>              | XtenCHO   |
| <b>REACTIVITY</b>        | Severe acute respiratory syndrome coronavirus 2 (2019-nCoV) (SARS-CoV-2)  |
| <b>FORM &amp; BUFFER</b> | 0.01M PBS, pH 7.4.  |
| <b>STORAGE</b>           | Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles. |

