



Proteogenix

Proteogenix SAS

www.proteogenix.science
19 rue de la Haye
67300 Schiltigheim, France

Tel +33 (0) 390 20 54 70
Fax +33 (0) 978 53 36 90
products@proteogenix.fr

Anti-SARS-CoV-2 S Protein Nanobody (SAA0996)

OVERVIEW

PRODUCT REFERENCE	PTX19025
CLONALITY	Monoclonal Antibody
ISOTYPE	VHH-8His-Cys-tag
TARGET	S Protein

PRODUCT INFORMATION

HOST	XtenCHO
REACTIVITY	Severe acute respiratory syndrome coronavirus 2 (2019-nCoV) (SARS-CoV-2)
FORM & BUFFER	0.01M PBS, pH 7.4.
STORAGE	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.



WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.