



## Anti-SARS-CoV-2 S Protein Nanobody (SAA1048)

### OVERVIEW

---

<b>PRODUCT REFERENCE</b>	PTX19030
<b>CLONALITY</b>	Monoclonal Antibody
<b>ISOTYPE</b>	VHH-8His-Cys-tag
<b>TARGET</b>	S Protein

### PRODUCT INFORMATION

---

<b>HOST</b>	XtenCHO
<b>REACTIVITY</b>	Severe acute respiratory syndrome coronavirus 2 (2019-nCoV) (SARS-CoV-2)
<b>FORM &amp; BUFFER</b>	0.01M PBS, pH 7.4.
<b>STORAGE</b>	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.

