



Anti-SARS-CoV-2 S Protein Nanobody (SAA1094)

OVERVIEW

PRODUCT REFERENCE	PTX19060
CLONALITY	Monoclonal Antibody
ISOTYPE	VHH-8His-Cys-tag
TARGET	S Protein

PRODUCT INFORMATION

HOST	XtenCHO
REACTIVITY	Severe acute respiratory syndrome coronavirus 2 (2019-nCoV) (SARS-CoV-2)
FORM & BUFFER	0.01M PBS, pH 7.4.
STORAGE	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.

