



Anti-SARS-CoV-2 S Protein Nanobody (SAA1102)

OVERVIEW

| | |
|--------------------------|---------------------|
| PRODUCT REFERENCE | PTX19065 |
| CLONALITY | Monoclonal Antibody |
| ISOTYPE | VHH-8His-Cys-tag |
| TARGET | S Protein |

PRODUCT INFORMATION

| | |
|--------------------------|---|
| HOST | XtenCHO |
| REACTIVITY | Severe acute respiratory syndrome coronavirus 2 (2019-nCoV) (SARS-CoV-2) |
| FORM & BUFFER | 0.01M PBS, pH 7.4. |
| STORAGE | Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles. |

