



# Proteogenix

**Proteogenix SAS**

www.proteogenix.science  
19 rue de la Haye  
67300 Schiltigheim, France  
Tel +33 (0) 390 20 54 70  
Fax +33 (0) 978 53 36 90  
products@proteogenix.fr

## Anti-HIV1 Surface protein gp120 Nanobody (SAA0912)

### OVERVIEW

---

|                          |                       |
|--------------------------|-----------------------|
| <b>PRODUCT REFERENCE</b> | PTX19079              |
| <b>CLONALITY</b>         | Monoclonal Antibody   |
| <b>ISOTYPE</b>           | VHH-8His-Cys-tag      |
| <b>TARGET</b>            | Surface protein gp120 |

### PRODUCT INFORMATION

---

|                          |   |
|--------------------------|---|
| <b>HOST</b>              | XtenCHO   |
| <b>REACTIVITY</b>        | Human immunodeficiency virus type 1 group M subtype B (isolate HXB2) (HIV-1)  |
| <b>FORM &amp; BUFFER</b> | 0.01M PBS, pH 7.4.  |
| <b>STORAGE</b>           | Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles. |



**WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.**