



Proteogenix

Proteogenix SAS

www.proteogenix.science
19 rue de la Haye
67300 Schiltigheim, France
Tel +33 (0) 390 20 54 70
Fax +33 (0) 978 53 36 90
products@proteogenix.fr

Anti-HIV1 Surface protein gp120 Nanobody (SAA0894)

OVERVIEW

PRODUCT REFERENCE	PTX19083
CLONALITY	Monoclonal Antibody
ISOTYPE	VHH-8His-Cys-tag
TARGET	Surface protein gp120

PRODUCT INFORMATION

HOST	XtenCHO
REACTIVITY	Human immunodeficiency virus type 1 group M subtype B (isolate HXB2) (HIV-1)
FORM & BUFFER	0.01M PBS, pH 7.4.
STORAGE	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.



WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.