



Proteogenix

Proteogenix SAS

www.proteogenix.science
19 rue de la Haye
67300 Schiltigheim, France
Tel +33 (0) 390 20 54 70
Fax +33 (0) 978 53 36 90
products@proteogenix.fr

Anti-IL1RL1/ST2 Monoclonal Antibody, C2244

OVERVIEW

PRODUCT REFERENCE	PTX19085
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1
TARGET	T1, DER4, Interleukin-1 receptor-like 1, ST2, IL1RL1, Protein ST2

PRODUCT INFORMATION

HOST	Human
REACTIVITY	Human
FORM & BUFFER	0.01M PBS, pH 7.4.
STORAGE	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.



WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.