



Anti-IL1RL1/ST2 Monoclonal Antibody, STLM208

OVERVIEW

PRODUCT REFERENCE	PTX19086
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1
TARGET	T1, DER4, Interleukin-1 receptor-like 1, ST2, IL1RL1, Protein ST2

PRODUCT INFORMATION

HOST	Human
REACTIVITY	Human
FORM & BUFFER	0.01M PBS, pH 7.4.
STORAGE	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.

