



Anti-Human Abeta42/Beta-APP42 Monoclonal Antibody, 21F12

OVERVIEW

PRODUCT REFERENCE	PTX19087
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1
TARGET	Abeta42, Beta-APP42, Amyloid-beta protein 42, amyloid beta A4 precursor protein, Amyloid-beta A4 protein

PRODUCT INFORMATION

HOST	Mouse
REACTIVITY	Human
FORM & BUFFER	0.01M PBS, pH 7.4.
STORAGE	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.

