



Proteogenix

Proteogenix SAS

www.proteogenix.science
19 rue de la Haye
67300 Schiltigheim, France

Tel +33 (0) 390 20 54 70
Fax +33 (0) 978 53 36 90
products@proteogenix.fr

Anti-DENV-2 E-DII/Envelope Protein E Monoclonal Antibody, 2A10

OVERVIEW

PRODUCT REFERENCE	PTX19089
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1, kappa
TARGET	Genome polyprotein, Envelope protein E, E protein domain II, E-DII, Dengue Virus (DENV)

PRODUCT INFORMATION

HOST	Human
REACTIVITY	Dengue Virus (DENV1-4), Yellow Fever Virus (YFV), Zika virus (ZIKV), West Nile virus (WNV)
FORM & BUFFER	0.01M PBS, pH 7.4.
STORAGE	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.



WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.