



## Anti-HLA-A1-peptide (EADPTGHSY) Monoclonal Antibody, HYB3

### OVERVIEW

---

<b>PRODUCT REFERENCE</b>	PTX19093
<b>CLONALITY</b>	Monoclonal Antibody
<b>ISOTYPE</b>	IgG1
<b>TARGET</b>	HLAA, Human leukocyte antigen A, HLA class I histocompatibility antigen, A alpha chain, HLA-A, MAGEA1, CT1.1, Cancer/testis antigen 1.1, MAGE1A, MAGE1, Antigen MZ2-E, Melanoma-associated antigen 1, MAGE-1 antigen

### PRODUCT INFORMATION

---

<b>HOST</b>	Mouse
<b>REACTIVITY</b>	Human
<b>FORM &amp; BUFFER</b>	0.01M PBS, pH 7.4.
<b>STORAGE</b>	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.

