



## Anti-HLA-A2-peptide (RMFPNAPYL) Monoclonal Antibody, ESK1

### OVERVIEW

---

<b>PRODUCT REFERENCE</b>	PTX19094
<b>CLONALITY</b>	Monoclonal Antibody
<b>ISOTYPE</b>	IgG1
<b>TARGET</b>	HLAA, Human leukocyte antigen A, HLA class I histocompatibility antigen, A alpha chain, HLA-A, WT1, Wilms tumor protein, WT33

### PRODUCT INFORMATION

---

<b>HOST</b>	Mouse
<b>REACTIVITY</b>	Human
<b>FORM &amp; BUFFER</b>	0.01M PBS, pH 7.4.
<b>STORAGE</b>	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.

