



## Anti-HLA-A2-peptide (LLFGYPVYV) Monoclonal Antibody, T3F2

### OVERVIEW

---

<b>PRODUCT REFERENCE</b>	PTX19096
<b>CLONALITY</b>	Monoclonal Antibody
<b>ISOTYPE</b>	IgG1
<b>TARGET</b>	HLAA, Human leukocyte antigen A, HLA class I histocompatibility antigen, A alpha chain, HLA-A, Protein Tax-1, Protein X-LOR, Protein PX, Trans-activating transcriptional regulatory protein of HTLV-1, Tax

### PRODUCT INFORMATION

---

<b>HOST</b>	Mouse
<b>REACTIVITY</b>	Human
<b>FORM &amp; BUFFER</b>	0.01M PBS, pH 7.4.
<b>STORAGE</b>	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.

