



## Anti-HLA-A2-peptide (ILAKFLHWL) Monoclonal Antibody, 4B4

### OVERVIEW

---

<b>PRODUCT REFERENCE</b>	PTX19098
<b>CLONALITY</b>	Monoclonal Antibody
<b>ISOTYPE</b>	IgG1
<b>TARGET</b>	HLAA, Human leukocyte antigen A, HLA class I histocompatibility antigen, A alpha chain, HLA-A, EST2, TCS1, HEST2, Telomerase reverse transcriptase, Telomerase catalytic subunit, Telomerase-associated protein 2, TRT, TP2, TERT

### PRODUCT INFORMATION

---

<b>HOST</b>	Mouse
<b>REACTIVITY</b>	Human
<b>FORM &amp; BUFFER</b>	0.01M PBS, pH 7.4.
<b>STORAGE</b>	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.

