



Anti-HLA-A2-peptide (ILAKFLHWL) Monoclonal Antibody, 4C2

OVERVIEW

PRODUCT REFERENCE	PTX19099
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1
TARGET	HLAA, Human leukocyte antigen A, HLA class I histocompatibility antigen, A alpha chain, HLA-A, EST2, TCS1, HEST2, Telomerase reverse transcriptase, Telomerase catalytic subunit, Telomerase-associated protein 2, TRT, TP2, TERT

PRODUCT INFORMATION

HOST	Mouse
REACTIVITY	Human
FORM & BUFFER	0.01M PBS, pH 7.4.
STORAGE	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.

