



Proteogenix

Proteogenix SAS

www.proteogenix.science
19 rue de la Haye
67300 Schiltigheim, France
Tel +33 (0) 390 20 54 70
Fax +33 (0) 978 53 36 90
products@proteogenix.fr

Anti-HLA-A2-peptide (RLVDDFLLV) Monoclonal Antibody, 3A12

OVERVIEW

PRODUCT REFERENCE	PTX19100
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1
TARGET	HLAA, Human leukocyte antigen A, HLA class I histocompatibility antigen, A alpha chain, HLA-A, EST2, TCS1, HEST2, Telomerase reverse transcriptase, Telomerase catalytic subunit, Telomerase-associated protein 2, TRT, TP2, TERT

PRODUCT INFORMATION

HOST	Mouse
REACTIVITY	Human
FORM & BUFFER	0.01M PBS, pH 7.4.
STORAGE	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.



WARNING: THIS PRODUCT IS FOR RESEARCH USE ONLY.