



Anti-LASV GPC Monoclonal Antibody, 37.2G

OVERVIEW

PRODUCT REFERENCE	PTX19106
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1, kappa
TARGET	Pre-glycoprotein polyprotein GP complex, Pre-GP-C, GPC, GP-C, Lassa Virus (LASV)

PRODUCT INFORMATION

HOST	Human
REACTIVITY	Lassa virus (strain Mouse/Sierra Leone/Josiah/1976) (LASV)
FORM & BUFFER	0.01M PBS, pH 7.4.
STORAGE	Store at 4 °C for short term only. Store between -20 °C and -80 °C for long term and thaw on ice (at 4 °C) before use. Avoid repeated freeze-thaw cycles.

